

ORDER INFORMATION ALVISION™ (continued)

		KAPLAN		KAPLAN JS	
		KAPLAN 4.0 PP	KAPLAN 4.5 PP	KAPLAN JS 3.5 PP	KAPLAN JS 4.0 PP
5 Pack	Length	100 cm	100 cm	100 cm	100 cm
	∅				
	5 F	501110019301	501110019302	501110019304	501110019303
5 Pack	Length	110 cm	110 cm	110 cm	110 cm
	∅				
	5 F	501110019309	501110019310	501110019312	501110019311
5 Pack	6 F	501110019313	501110019314	501110019316	501110019315

		ATESAL		
		ATESAL 3.5	ATESAL 4.0	ATESAL 4.5
5 Pack	Length	100 cm	100 cm	100 cm
	∅			
	4 F	501110019641	501110019640	501110019642
5 Pack	Length	110 cm	110 cm	110 cm
	∅			
	5 F	501110019740	501110019741	501110019742
5 Pack	6 F	501110019840	501110019841	501110019842
5 Pack	7 F*	501110019941	501110019940	501110019942

		KAPLAN MOD I		
		KAPLAN MOD I 4.0 PP	KAPLAN MOD I 4.5 PP	KAPLAN MOD I 5.0 PP
5 Pack	Length	100 cm	100 cm	100 cm
	∅			
	5 F	501110019786	501110019787	501110019890
5 Pack	Length	110 cm	110 cm	110 cm
	∅			
	5 F	501110019893	501110019894	501110019892
5 Pack	6 F	501110019793	501110019794	501110019792

		KAPLAN MOD II		
		KAPLAN MOD II 4.0 PP	KAPLAN MOD II 4.5 PP	KAPLAN MOD II 5.0 PP
5 Pack	Length	100 cm	100 cm	100 cm
	∅			
	5 F	501110019886	501110019887	501110019891
5 Pack	Length	110 cm	110 cm	110 cm
	∅			
	5 F	501110019896	501110019897	501110019895
5 Pack	6 F	501110019796	501110019797	501110019795

* Check with your sales representative the availability

PP = PRO PERFUSION - CURVES WITH SIDEHOLES

Kaplan MOD I and II : one frontal side hole
Kaplan PP: one frontal side hole
Kaplan JS PP: two lateral side holes



2023/06 1C70155095/4

ORDER INFORMATION ALVISION™

		CORONARY BYPASS		CORONARY BYPASS GRAFT		INTERNAL MAMMARY	
		CORONARY BYPASS LEFT (LCB)	CORONARY BYPASS RIGHT (RCB)	CORONARY BYPASS GRAFT LEFT	CORONARY BYPASS GRAFT RIGHT	INTERNAL MAMMARY (IMA)	
5 Pack	Length	100 cm	100 cm	100 cm	100 cm	100 cm	
	∅						
	4 F	501110019621	501110019622	501110019652	501110019651	501110019623	
5 Pack	Length	110 cm	110 cm	110 cm	110 cm	110 cm	
	∅						
	5 F	501110019721	501110019722	501110019752	501110019751	501110019723	
5 Pack	Length	110 cm	110 cm	110 cm	110 cm	110 cm	
	∅						
	6 F	501110019821	501110019822	501110019852	501110019851	501110019823	
5 Pack	7 F*	501110019921	501110019922	501110019952	501110019951	501110019923	

* Check with your sales representative the availability

TECHNICAL DATA

Category:	Interventional cardiology diagnostic catheters
Diameters:	4F, 5F, 6F and 7F
Shaft :	Multi-segments (with variable stiffness)
Proximal shaft material:	Braided Nylon
Braiding:	Stainless steel flat wires
Distal shaft material:	Radiopaque PEBA
Tip material:	Soft PEBA with radiopaque Barium Sulfate
Guidewire compatibility:	0.038"
Maximum pressure:	1200psi

Flow Rate

Size	Catheter			Pressure			
	Curve	ID	Length	300psi	600psi	900psi	1200psi
4F	Selective	0.042"	100 cm	11ml/s	15ml/s	18ml/s	22ml/s
	Pigtail	0.042"	110 cm	11ml/s	15ml/s	19ml/s	22ml/s
	Selective	0.042"	125 cm	10ml/s	14ml/s	17ml/s	20ml/s
	Pigtail	0.042"	125 cm	10ml/s	14ml/s	17ml/s	21ml/s
5F	Selective	0.042"	80 cm	12ml/s	16ml/s	20ml/s	24ml/s
	Selective	0.047"	100 cm	14ml/s	19ml/s	24ml/s	28ml/s
	Pigtail	0.047"	110 cm	14ml/s	19ml/s	25ml/s	28ml/s
	Selective	0.047"	110 cm	13ml/s	18ml/s	24ml/s	26ml/s
6F	Selective	0.047"	125 cm	13ml/s	17ml/s	22ml/s	25ml/s
	Pigtail	0.047"	125 cm	13ml/s	18ml/s	23ml/s	26ml/s
	Selective	0.047"	80 cm	15ml/s	21ml/s	26ml/s	30ml/s
	Selective	0.057"	100 cm	23ml/s	33ml/s	39ml/s	41ml/s
7F	Pigtail	0.054"	110 cm	21ml/s	29ml/s	35ml/s	37ml/s
	Selective	0.057"	110 cm	23ml/s	32ml/s	38ml/s	41ml/s
	Selective	0.057"	125 cm	21ml/s	30ml/s	35ml/s	37ml/s
	Pigtail	0.054"	125 cm	19ml/s	27ml/s	33ml/s	36ml/s
7F	Selective	0.057"	80 cm	26ml/s	36ml/s	41ml/s	45ml/s
	Selective	0.070"	100 cm	38ml/s	45ml/s	45ml/s	45ml/s
	Pigtail	0.067"	110 cm	35ml/s	45ml/s	45ml/s	45ml/s
	Selective	0.070"	80 cm	43ml/s	45ml/s	45ml/s	45ml/s

www.alvimedica.com

Alvision™ Interventional Cardiology Diagnostic Catheter

Please contact your local sales representative for product availability.
The images shown within this document are of a purely illustrative nature and are not intended to depict reality.

Please consult product labels and package inserts for indications, contraindications, hazards, warnings, cautions and instruction for use.
©2023, Alvimedica. All Rights Reserved.



Alvision™


INTERVENTIONAL CARDIOLOGY DIAGNOSTIC CATHETER


Coronary angiographic catheters designed for effective diagnostic procedures

- Outstanding torque response and catheter stability due to Flat Wire Braid design
- Large inner lumens allow for superior visualization
- Wide choice of shapes to master every challenge



ORDER INFORMATION ALVISION™ (continued)


JUDKINS LEFT							
							
	JUDKINS LEFT (JL 3.0)	JUDKINS LEFT (JL 3.5)	JUDKINS LEFT (JL 4.0)	JUDKINS LEFT (JL 4.5)	JUDKINS LEFT (JL 5.0)	JUDKINS LEFT (JL 6.0)	JUDKINS LEFT ST (JL 4.0 ST)
Length	100 cm	100 cm	100 cm	100 cm	100 cm	100 cm	100 cm
Ø							
5 Pack	4 F 501110019601	501110019602	501110019603	501110019605	501110019606	501110019607	501110019604
	5 F 501110019701	501110019702	501110019704	501110019705	501110019706	501110019707	501110019703
	6 F 501110019801	501110019802	501110019803	501110019805	501110019806	501110019807	501110019804
	7 F* 501110019901	501110019902	501110019903	501110019905	501110019906	501110019907	501110019904
Length			125 cm		125 cm	125 cm	
Ø							
5 Pack	4 F		501110019502		501110019503	501110019501	
	5 F		501110019507		501110019508	501110019509	
	6 F		501110019513		501110019514	501110019515	


JUDKINS RIGHT						
						
	JUDKINS RIGHT (JR 3.0)	JUDKINS RIGHT (JR 3.5)	JUDKINS RIGHT (JR 4.0)	JUDKINS RIGHT (JR 4.5)	JUDKINS RIGHT (JR 5.0)	JUDKINS RIGHT (JR 6.0)
Length	100 cm	100 cm	100 cm	100 cm	100 cm	100 cm
Ø						
5 Pack	4 F 501110019608	501110019609	501110019610	501110019611	501110019612	501110019613
	5 F 501110019708	501110019709	501110019711	501110019710	501110019712	501110019713
	6 F 501110019808	501110019809	501110019810	501110019811	501110019812	501110019813
	7 F* 501110019908	501110019909	501110019910	501110019911	501110019912	501110019913
Length			125 cm		125 cm	125 cm
Ø						
5 Pack	4 F		501110019504		501110019505	501110019506
	5 F		501110019510		501110019511	501110019512
	6 F		501110019516		501110019517	501110019518


* Check with your sales representative the availability

PP = PRO PERFUSION - CURVES WITH SIDEHOLES ST = SHORT TIP

ORDER INFORMATION ALVISION™ (continued)

AMPLATZ LEFT			
			
	AMPLATZ LEFT (AL 1)	AMPLATZ LEFT (AL 2)	AMPLATZ LEFT (AL 3)
Length	100 cm	100 cm	100 cm
Ø			
5 Pack	4 F 501110019614	501110019615	501110019616
	5 F 501110019714	501110019715	501110019716
	6 F 501110019814	501110019815	501110019816
	7 F* 501110019914	501110019915	501110019916


AMPLATZ RIGHT				
				
	AMPLATZ RIGHT (AR 1)	AMPLATZ RIGHT (AR 2)	AMPLATZ RIGHT (AR 3)	AMPLATZ RIGHT MODIFIED B (AR MOD B)
Length	100 cm	100 cm	100 cm	100 cm
Ø				
5 Pack	4 F 501110019617	501110019618	501110019619	501110019620
	5 F 501110019717	501110019718	501110019719	501110019720
	6 F 501110019817	501110019818	501110019819	501110019820
	7 F* 501110019917	501110019918	501110019919	501110019920





PIGTAIL			
			
	PIGTAIL STRAIGHT (PIG)	PIGTAIL 145 (PIG 145)	PIGTAIL 155 (PIG 155)
Length	110 cm	110 cm	110 cm
Ø			
5 Pack	4 F 501110019637	501110019638	501110019639
	5 F 501110019737	501110019738	501110019739
	6 F 501110019837	501110019838	501110019839
	7 F* 501110019937	501110019938	501110019939
Length	125 cm	125 cm	125 cm
Ø			
5 Pack	4 F 501110019525	501110019519	501110019520
	5 F 501110019526	501110019521	501110019522
	6 F 501110019527	501110019523	501110019524

* Check with your sales representative the availability

PP = PRO PERFUSION - CURVES WITH SIDEHOLES ST = SHORT TIP

ORDER INFORMATION ALVISION™ (continued)

SONES						
						
	SONES I (SON 1)	SONES II (SON 2)	SONES III (SON 3)	SONES I (SON 1)	SONES II (SON 2)	SONES III (SON 3)
Length	100 cm	100 cm	100 cm	80 cm	80 cm	80 cm
Ø						
5 Pack	4 F 501110019628	501110019629	501110019630	501110019631	501110019632	501110019633
	5 F 501110019728	501110019729	501110019730	501110019731	501110019732	501110019733
	6 F 501110019828	501110019829	501110019830	501110019831	501110019832	501110019833
	7 F* 501110019928	501110019929	501110019930	501110019931	501110019932	501110019933
Length						
Ø						
5 Pack	PP	SONES I PP (SON 1 PP)	SONES II PP (SON 2 PP)	SONES III PP (SON 3 PP)		
	4 F 501110019634	501110019635	501110019636			
	5 F 501110019734	501110019735	501110019736			
	6 F 501110019834	501110019835	501110019836			
	7 F* 501110019934	501110019935	501110019936			

Q4		HOCKEY STICK		MOD. EXTRA BACK UP		PROGRESSIVE RIGHT	
							
	Q4		HOCKEY STICK		MOD. EXTRA BACK UP		PROGRESSIVE RIGHT
Length	100 cm		100 cm		100 cm		100 cm
Ø							
5 Pack	4 F 501110019647		4 F 501110019643		4 F 501110019648		4 F 501110019644
	5 F 501110019747		5 F 501110019743		5 F 501110019748		5 F 501110019744
	6 F 501110019847		6 F 501110019843		6 F 501110019848		6 F 501110019844
	7 F* 501110019947		7 F* 501110019943		7 F* 501110019948		7 F* 501110019944

* Check with your sales representative the availability

PP = PRO PERFUSION - CURVES WITH SIDEHOLES ST = SHORT TIP